


QuickFISH® Fixation Kit - BC

QuickFISH Fixation Kit for Blood Cultures



IVD

REF

QFFIXBC1-25

Contents

The reagents and slides included in the QuickFISH BC Fixation Kit are for use with the QuickFISH BC culture identification kits according to the instructions outlined in the respective package inserts.

Reagents

QuickFISH BC Fixation Kit

Each kit contains sufficient material for 25 tests. Reagents are supplied ready for use. The expiration date of the kit is as indicated on the outer box label.

QuickFISH Slide QuickFISH Slide CS012

25 QuickFISH slides with controls.

QuickFix-1 QuickFix-1 CP0169




Primary fixation solution.

QuickFix-2 QuickFix-2 CP0170

Secondary fixation solution.

Warnings and Precautions

Safety Precautions

QuickFix-1	Contains 24% ethanol	Harmful to aquatic life with long lasting effects. Safety Data Sheet is available upon request. Available in the QuickFISH Fixation Kit.
QuickFix-2	   Danger Contains 97% methanol	Highly flammable liquid and vapor. Toxic if swallowed. Tonic in contact with skin. Toxic if inhaled. Causes damage to the central nervous system. Safety Data Sheet is available upon request. Available in the QuickFISH Fixation Kit.

Do not eat, drink, smoke, apply cosmetics, store or prepare foods within the designated work area.

Dispose of reagents in accordance with federal, state, and local regulations.

Technical Precautions

Reagents must not be used after the expiration dates printed on the labels.

Reagents are provided at fixed concentrations. Assay performance may be affected if the reagents are modified in any way or are not stored

under the recommended conditions as detailed in "Storage of Kit Components".

Avoid microbial contamination of reagents.

Avoid any cross-contamination of samples and reagents, as this may give rise to erroneous results.

Do not allow dropper bottle tip to touch the smear as this may cause cross contamination of material between slides, or cause contamination of the reagent.









Storage and Preparation of Kit Components

To ensure optimal kit performance, it is important that kit components are stored according to the following instructions:

Store kit components at 2-8°C. Store bottles upright and tighten caps after use. Place kit components at room temperature prior to use and return the kit components to 2-8°C after use.

QuickFISH slides are provided in individually sealed pouches with nitrogen and a desiccant. Store slides at 2-8°C. Slides must be used immediately after breaking pouch seal. Do not use slides after the expiration date.

Definitions

REF	Product code/catalog number	LOT	Batch code
	Consult the instructions for use		Storage temperature limitations
	Contains sufficient for <n> tests		Health Hazard
	Manufacturer		Skull and Crossbones
EC REP	Authorized Representative		Flame
	Use by		

Technical Advice and Customer Service

For all inquiries, please contact OpGen or your local distributor.



OpGen, Inc.
708 Quince Orchard Rd
Gaithersburg, MD 20878
USA

Tel: +1 301 869 9683
Fax: +1 301 869 9684

techsupport@opgen.com



Curetis GmbH
Max-Eyth-Straße 42
71088 Holzgerlingen,
Germany

Tel: +49 7031 49195 10
Fax: +49 7031 49195 19

www.OpGen.com

30 April 2020

PN2048E-EN
DCR 20-0035