



## Unyvero LRT/Pneumonia Control Panel M404

The M404 kit includes 12 single-use vials of ready-to-use qualitative external controls with coverage for all targets on the Unyvero LRT Cartridge. Each kit contains:

- 4 tubes x Unyvero LRT/Pneumonia Positive A
- 4 tubes x Unyvero LRT/Pneumonia Positive B
- 4 tubes x Unyvero LRT/Pneumonia Negative

**MMQCI Part Number:** M404

**Storage Conditions (M404 & M408):** Refrigerated 2 – 8°C. Unopened tubes are stable through the expiration date printed on the kit label when appropriately stored.

## Unyvero LRT Verification Panel M408

The M408 kit includes 27 single-use vials of ready-to-use qualitative verification reference material with coverage for all targets on the Unyvero LRT Cartridge. Each kit contains:

- 4 tubes x Unyvero LRT Verification 1
- 4 tubes x Unyvero LRT Verification 2
- 4 tubes x Unyvero LRT Verification 3
- 4 tubes x Unyvero LRT Verification 4
- 4 tubes x Unyvero LRT Verification 5
- 7 tubes x Unyvero LRT Verification Negative

**MMQCI Part Number:** M408

### Ordering Information:

Contact MMQCI for all questions and orders. Submit orders by phone, fax or email:

phone: **1-207-885-1072 ext. 1**  
 fax: 1-207-885-1079  
 email: **info@mmqci.com**

